

Venus Versa User Manual (International Version)

Document part number: DCVEROW77B

October 2017

Copyright © 2017 Venus Concept Ltd. All rights reserved.

Venus Concept Ltd. reserves the right to make changes to its products or specifications to improve performance, reliability, or manufacturability. Information furnished by Venus Concept Ltd. is believed to be accurate and reliable. However, Venus Concept Ltd. assumes no responsibility for its use. No license is granted by its implication or otherwise under any patent or patent rights of Venus Concept Ltd.

No part of this document may be produced or transmitted in any form or by any means, electronic or mechanical, for any purpose, without a written permission of Venus Concept Ltd.

Data and Specifications are subject to change without notification.

Venus Concept Ltd. has patents and pending patent applications, trademarks, copyrights, or other intellectual property rights covering subject matter in this document. The furnishing of this document does not give you any license to these patents, trademarks, copyrights, or other intellectual property rights.

NOTICE: This document may contain Venus Concept Confidential Information. If you received this document by mistake, please inform us, so as to enable us to retrieve the documents. In such cases, copying and distributing this document are strictly prohibited.

Table of Contents

CHAPTER 1	About this Manual	Page
1.1.	Use of the Manual.....	1-2
1.2.	Scope of the Manual	1-3
1.3.	Document Conventions.....	1-6
1.4.	Symbols	1-7

CHAPTER 2	Safety and Regulatory Information	Page
2.1.	Introduction.....	2-2
2.2.	Intended Use	2-3
2.3.	Intended User and Environment.....	2-3
2.4.	Safety Information	2-4
2.5.	Electrical and Mechanical Safety	2-7
2.6.	Fire Safety.....	2-8
2.7.	Protective Eyewear	2-8
2.8.	Optical Safety during the IPL Treatment..	2-9
2.9.	RF Treatment Safety Precautions	2-10
2.10.	Compliance with International Standards ...	2-11
2.11.	System Safety Features.....	2-12
2.12.	Electromagnetic Compatibility	2-13
2.13.	System Labels.....	2-18
2.14.	Contact Information.....	2-20

CHAPTER 3	Installation	Page
3.1.	Installing the System	3-2
3.2.	Electrical Requirements	3-2
3.3.	Environmental Requirements	3-3
3.4.	Equipment List	3-4
3.5.	Connecting the System to the Mains Power.....	3-5
3.6.	Connecting the Applicators	3-6
3.6.1.	Connecting the IPL Applicator.....	3-6
3.6.2.	Connecting the RF Applicators	3-7
3.6.3.	Connecting the Viva Applicator Tip.....	3-8

CHAPTER 4	System Description	Page
4.1.	General Description.....	4-2
4.2.	System Components	4-3
4.2.1.	System Console	4-4
4.2.2.	Touch Screen	4-4
4.2.3.	Emergency Button	4-5
4.2.4.	IPL Applicator Connection Receptacle.....	4-6
4.2.5.	RF Applicator Connection Receptacles.....	4-7
4.2.6.	IPL Applicator Holder.....	4-8
4.2.7.	RF Applicator Holders	4-9
4.2.8.	Main Power Switch.....	4-10
4.2.9.	Power Connection Port.....	4-11

CHAPTER 4	System Description (cont.)	Page
4.3.	IPL Applicators.....	4-12
4.3.1.	SR515 Applicator.....	4-12
4.3.2.	SR580 Applicator.....	4-13
4.3.3.	HR650 Applicator.....	4-14
4.3.4.	HR690 Applicator.....	4-15
4.3.5.	HR650XL Applicator	4-16
4.3.6.	HR690XL Applicator	4-17
4.3.7.	ACDUAL Applicator.....	4-18
4.4.	RF Applicators.....	4-19
4.4.1.	Diamondpolar Applicator.....	4-19
4.4.2.	Octipolar Applicator.....	4-20
4.4.3.	Viva Applicator.....	4-21

CHAPTER 5	Operating Instructions	Page
5.1.	Turning the System On/Off	5-2
5.2.	Logging In.....	5-4
5.3.	Selecting the Desired Applicator	5-6
5.4.	Treatment Modes and Parameters	5-7
5.4.1.	(MP) ² Treatment Mode	5-8
5.4.2.	Fractional RF Treatment Mode	5-12
5.4.3.	IPL Treatment Mode.....	5-17
5.5.	Setting the Treatment Parameters.....	5-21
5.6.	Using Preset Parameters in the (MP) ² Treatment Mode.....	5-22

CHAPTER 5	Operating Instructions (cont.)	Page
	5.6.1. <i>Recalling the Preset Treatment Parameters</i>	5-22
	5.6.2. <i>Changing the Preset Treatment Parameters</i>	5-23
5.7.	Emitting the Treatment Energy	5-25

CHAPTER 6	Advanced Features and Settings	Page
6.1.	Tools screen.....	6-2
	6.1.1. <i>Features in the Tools Screen</i>	6-2
	6.1.2. <i>Status Bar</i>	6-5
	6.1.3. <i>Accessing the Tools Screen</i>	6-5
6.2.	Setting the Interface language	6-6
6.3.	Setting the Volume	6-7
6.4.	Setting the Login Code	6-8
6.5.	Logging Out.....	6-11
6.6.	Setting the Auto Logout Time	6-12
6.7.	Displaying the System Information	6-13
6.8.	Extending the Activation Period	6-14
6.9.	Displaying the System Counter Information	6-16
6.10.	Activating the Pump	6-18
6.11.	Technician Mode	6-18

CHAPTER 7	Care and Maintenance	Page
7.1.	Cleaning the System	7-2
7.2.	Cleaning the Applicators	7-2
7.3.	Disinfecting the Applicators	7-3
7.4.	Disinfecting the Viva Applicator's Tip	7-3
7.5.	Sterilizing the Viva Applicator's Tip	7-4
7.6.	Adding Water to the Cooling System	7-5
7.7.	Circulating the Water	7-10
7.8.	Replacing a Burnt Fuse	7-15
7.9.	Moving the System	7-16

CHAPTER 8	Troubleshooting	Page
8.1.	General	8-2
8.2.	System Problems	8-3
8.3.	Errors and Messages	8-5

CHAPTER 9	Specifications	Page
9.1.	System Specifications.....	9-2
9.2.	SR515 Applicator Specifications	9-3
9.3.	SR580 Applicator Specifications	9-4
9.4.	HR650 Applicator Specifications.....	9-5
9.5.	HR690 Applicator Specifications.....	9-6
9.6.	HR650XL Applicator Specifications	9-7
9.7.	HR690XL Applicator Specifications	9-8
9.8.	ACDUAL Applicator Specifications.....	9-9
9.9.	Diamondpolar Applicator Specifications	9-10
9.10.	Octipolar Applicator Specifications	9-11
9.11.	Viva Applicator Specifications	9-12
9.12.	Environmental Requirements	9-13
9.13.	Output Power Graphs	9-14

Appendix A	Clinical Guide for (MP) ² Treatment with Diamondpolar Applicator	Page
A.1.	Introduction.....	A-2
A.2.	Indications.....	A-2
A.3.	Contraindications.....	A-3
A.4.	Exclusion Criteria	A-3
A.5.	Possible Side Effects	A-4
A.6.	Precaution Criteria	A-4
A.7.	Pre-Treatment Information	A-5
	<i>A.7.1. Recommended Accessories</i>	<i>A-5</i>
	<i>A.7.2. Clinical Data Management.....</i>	<i>A-5</i>
	<i>A.7.3. Expected Clinical Effects.....</i>	<i>A-6</i>
	<i>A.7.4. Adverse Effect</i>	<i>A-6</i>
	<i>A.7.5. Documentation.....</i>	<i>A-6</i>
	<i>A.7.6. Preparing for the Treatment.....</i>	<i>A-7</i>
A.8.	General Treatment Information	A-7
A.9.	Post-Treatment Care	A-9

Appendix B	Clinical Guide for (MP) ² Treatment with Octipolar Applicator	Page
B.1.	Introduction	B-2
B.2.	Indications	B-2
B.3.	Exclusion Criteria.....	B-3
B.4.	Precaution Criteria.....	B-4
B.5.	Pre-Treatment Information.....	B-5
	<i>B.5.1. Recommended Accessories.....</i>	<i>B-5</i>
	<i>B.5.2. Clinical Data Management</i>	<i>B-5</i>
	<i>B.5.3. Expected Clinical Effects.....</i>	<i>B-6</i>
	<i>B.5.4. Adverse Effect.....</i>	<i>B-6</i>
	<i>B.5.5. Documentation</i>	<i>B-7</i>
	<i>B.5.6. Preparing for the Treatment.....</i>	<i>B-8</i>
B.6.	General Treatment Information.....	B-8

Appendix C	Clinical Guide for Fractional RF Treatment with Viva Applicator	Page
C.1.	Introduction.....	C-2
C.2.	Indications.....	C-2
C.3.	Contraindications.....	C-3
C.4.	Exclusion Criteria	C-3
C.5.	Possible Side Effects	C-6
C.6.	Pre-Treatment Information	C-7
C.7.	Preparing the Applicator and Tip	C-8
C.8.	Test Treatment.....	C-9
C.9.	Treatment Procedure.....	C-10
C.10.	Treatment Protocol	C-13
C.11.	Treatment Parameters	C-13
C.12.	Post Treatment Care	C-14
C.13.	Follow-Up.....	C-15
C.14.	Treatment Conclusion.....	C-15

Appendix D	Clinical Guide for SR515 and SR580 Applicators	Page
D.1.	Introduction	D-2
D.2.	Indications	D-3
D.3.	Contraindications	D-4
D.4.	Precautions	D-5
D.5.	Possible Side Effects	D-5
D.6.	Before Starting	D-6
D.7.	Test Spot	D-7
D.8.	Treatment Method	D-7
D.9.	Treatment Protocol	D-9
	<i>D.9.1. Setting the Treatment Parameters</i>	<i>D-9</i>
	D.9.1.1. Pigmented Lesions	D-10
	D.9.1.2. Vascular Lesions	D-11
D.10.	Post-Treatment Care	D-12

Appendix E	Clinical Guide for Hair Removal Using IPL HR Applicators	Page
E.1.	Introduction	E-2
E.2.	Indications	E-3
E.3.	Contraindications	E-4
E.4.	Precautions	E-5
E.5.	Possible Side Effects	E-5
E.6.	Before Starting	E-6
E.7.	Test Spot	E-7
E.8.	Treatment Method	E-7
E.9.	Treatment Protocol	E-9
	<i>E.9.1. Setting the Treatment Parameters</i>	<i>E-10</i>
E.10.	Post-Treatment Care	E-11

pendix F	Clinical Guide for the ACDUAL IPL Applicator	Page
F.1.	Introduction	F-2
F.2.	Indications	F-3
F.3.	Contraindications	F-4
F.4.	Precautions	F-5
F.5.	Possible Side Effects	F-5
F.6.	Before Starting	F-6
F.7.	Test Spot.....	F-7
F.8.	Treatment Method.....	F-7
F.9.	Treatment Protocol.....	F-9
F.9.1.	Setting the Treatment Parameters.....	F-9
F.10.	Post-Treatment Care.....	F-10

pendix G	Manufacturer Warranty	Page
G.1.	Manufacturer Warranty	G-2

List of Figures

Figure	Description	Page
2-1.	IPL Treatment Room Caution Sign	2-9
2-2.	System Identification Label.....	2-18
2-3.	Light Emission Warning Label	2-19
3-1.	Power Connection Port on the System's Rear Side	3-5
3-2.	Connecting the IPL Applicator's Cable Connector	3-6
3-3.	Connecting the RF Applicator's Cable Connector	3-7
3-4.	Connecting a Tip to the Viva Applicator	3-8
4-1.	System Components	4-3
4-2.	Emergency Button	4-5
4-3.	IPL Applicator Connection Receptacle.....	4-6
4-4.	RF Applicator Connection Receptacles	4-7
4-5.	IPL Applicator Holder.....	4-8
4-6.	RF Applicator Holders	4-9
4-7.	Main Power Switch	4-10
4-8.	Power Connection Port.....	4-11
4-9.	SR515 Applicator	4-12
4-10.	SR580 Applicator	4-13
4-11.	HR650 Applicator	4-14
4-12.	HR690 Applicator	4-15
4-13.	HR650XL Applicator	4-16
4-14.	HR690XL Applicator	4-17
4-15.	ACDUAL Applicator	4-18
4-16.	Diamondpolar Applicator.....	4-19
4-17.	Octipolar Applicator.....	4-20
4-18.	Viva Applicator	4-21
5-1.	Power Connection Port.....	5-2
5-2.	Login Screen.....	5-3
5-3.	Entering the Login Code	5-4

Figure	Description	Page
4.	Treatment Screen (with the Diamondpolar Applicator Selected)	5-5
5.	Applicator Selector (with the Diamondpolar Applicator Selected)	5-6
6.	Elements of the (MP) ² Mode Treatment Screen (with the Octipolar Applicator Selected)	5-8
7.	Elements of the Fractional RF Mode Treatment Screen	5-12
8.	Elements of the Fractional RF Mode Advanced Treatment Screen	5-12
9.	Elements of the IPL Mode Treatment Screen (with the HR690 IPL Applicator Selected)	5-17
5-10.	Preset Screen	5-22
5-11.	Preset Icon Inward Arrow Changing to Green	5-24
6-1.	Tools Screen	6-2
6-2.	Status Bar	6-5
6-3.	Volume Screen	6-7
6-4.	Login Screen	6-8
6-5.	Entering the New Login Code	6-9
6-6.	Re-Entering the New Login Code	6-10
6-7.	Login Screen	6-11
6-8.	Set Auto Logout Time Screen	6-12
6-9.	Information Screen	6-13
6-10.	Activation Code Screen	6-15
6-11.	System Counter Information Screen	6-17
7-1.	Water Level Indicator on the System's rear Side	7-5
7-2.	Water Filling Kit	7-6
7-3.	Disconnecting the IPL Applicator's Cable Connector	7-6
7-4.	Connecting Two Parts of the Water Filling Kit'	7-7
7-5.	Inserting the Water Filling Kit's Connector	7-7
7-6.	Locking the Water Filling Kit's Connector	7-8
7-7.	Water Filling Kit Locked in Place	7-8

Figure	Description	Page
7-8.	Water Filling Kit Connected to the System	7-9
7-9.	Bypass Tool (Without the Funnel)	7-10
7-10.	Bypass Tool Connected to the System (Without the Funnel)	7-11
7-11.	Tools Screen	7-12
7-12.	Connecting the Bypass Tool and Starting the Pump Message	7-12
7-13.	Pump Activation Screen	7-13
7-14.	Water Circulation Completion Message	7-14
7-15.	Power Connection Port and Main Fuses on the System's Rear Side	7-15
9-1.	Graph of Output Power vs. Load	9-14
9-2.	Graph of Output Power vs. Setting	9-14
A-1.	Diamondpolar Applicator	A-2
B-1.	Octipolar Applicator	B-2
C-1.	Viva Applicator	C-2
D-1.	SR515 (Left) and SR580 (Right) IPL Applicators	D-2
E-1.	HR650 (Upper Left), HR690 (UpperRight), HR650XL (Lower Left) and HR690XL (Lower Right) PL Applicators	E-2
F-1.	ACDUAL Applicator	F-2

List of Tables

Table	Description	Page
1-1.	Glossary of Symbols Used in this Manual or on the Product	1-7
2-1.	Manufacturer's declaration – Electromagnetic Emissions	2-13
2-2.	Manufacturer's Declaration – Electromagnetic Immunity	2-14
2-3.	Manufacturer's Declaration – Electromagnetic Immunity that is not Life-Supporting	2-15
2-4.	Recommended Separation Distances between Portable and Mobile RF Communications Equipment and the ME Equipment or ME System that is not Life-Supporting.....	2-17
6-1.	Tools Screen Features	6-3
8-1.	Troubleshooting Problems that are Not Associated with a Message	8-3
8-2.	Troubleshooting Problems that are Associated with a Message	8-5
D-1.	Treatment Parameters for Pigmented Lesions	D-10
D-2.	Treatment Parameters for Vascular Lesions	D-11
E-1.	Treatment Parameters for Hair Removal	E-10
F-1.	Treatment Parameters for Acne	F-9

Chapter 1

About this Manual

Chapter Contents:

Section	Title	Page
1.1.	Use of the Manual.....	1-2
1.2.	Scope of the Manual	1-3
1.3.	Document Conventions	1-6
1.4.	Symbols	1-7